

PACKAGING / SHIPPING INSTRUCTIONS

Specimen Transport Materials

-
- **Specimens for all test panels are shipped at room temperature.**
 - **Each biohazard bag may contain only one patient sample.**
-


1 Place the labeled transport tube in the zippered pouch in the front of the biohazard bag with the absorbent pad. Seal the bag.



2 Place the completed Diatherix Eurofins Test Requisition in the flap pocket (no zipper) in the back of the biohazard transport bag.

Verify that the sample tube and the requisition are in the correct pouch:

- Requisition is located in the back side in the flapped pouch
- Specimen sample is located in the front side in the zippered pouch



3 Place the biohazard specimen bag(s), each containing only one sample, in the cardboard shipping box.

More than one specimen bag may be shipped in the same box.



4 Place the cardboard shipping box into the FedEx® UN 3373 Clinical Pak and seal it.

More than one box may be placed in the Clinical Pak.



5 Apply the barcoded shipping label to the FedEx® UN 3373 Clinical Pak. Retain the confirmation label for your records.

6 Call **1-800-463-3339** (1-800-GOFEDEX) to schedule a pickup. Listen to the automated system until after it states “*in a few words please describe what you are calling about.*” Then, press “0” for the next three prompts to be connected to an operator. State that you have an *Express Prepaid Label* and provide your address. No other information will be needed.

If your package cannot be picked up, it may be placed in a FedEx® Drop Box. A FedEx® Drop Box can be located by calling FedEx® or going to FedEx.com. Do not take this package to a FedEx® Office location as they will not accept this type of shipment.



TIP: For your records, write the pickup confirmation number in the space provided on the confirmation label.